



# Initial™ LiSi Block and G-CEM ONE™

## The perfect pairing.

Manufacturing of dental restorations has become a lot more potent with CAD/CAM. And it will increase to do so with GC's solutions. Initial LiSi Block and G-CEM ONE bring speed and efficiency into your practice, without any compromise on strength or aesthetics.



The patient was experiencing discomfort during chewing and hypersensitivity to cold in the first quadrant.



Teeth 15,16 and 17 after removal of old restorations and caries.



The teeth were built up again with composite and prepped for minimally invasive lithium disilicate restorations.



With Initial LiSi Block, beautiful and strong restorations can be milled quickly.



The optional G-CEM ONE Adhesive Enhancing Primer accelerates the polymerization and enhances the self-curing ability.



G-CEM ONE is applied to the inner surface of the crown that was pretreated in accordance with the IFU.



Excess cement is cleaned up fast.



Naturally beautiful.

Courtesy of Dr Claudia Mazzitelli, Italy



GC EUROPE N.V.

info.gce@gc.dental

<https://europe.gc.dental/en-GB>