

Washable Model resin offers high precision and strength, making it an ideal material for various applications. Unlike other resins, it can be easily cleaned with water instead of alcohol, resulting in smooth and non-sticky surfaces after curing. This makes it an exceptional choice for orthodontic use.

Applications

- Full Arch Model
- Diagnostic Models
- Models for Aligners
- Orthodontics Models

Advantages



Water washable



No IPA



Rapid curing



Color:



ROI

32 printed models of a 1 kg bottle

23 g/set

Property	Procedure	Specification
Resin viscosity (cps @25°C)	ASTM D7867	100-350
Tensile strength (Mpa)	ASTM D638	38.7
Flexural strength (Mpa)	ASTM D790	56.4
Flexural modulus (Mpa)	ASTM D790	1683



kexcelled
Printing excellence

MOLEGRID™

Dental 3D Printing Materials

Reach out to our sales team: dental@kexcelled3d.com



Follow us on LinkedIn



Visit our website